

HiTrap Fibro PrismA units

HiScreen Fibro PrismA units

AFFINITY CHROMATOGRAPHY

More targeted patient populations and smaller manufacturing batch sizes are driving the demand for increased process efficiency and flexible multiproduct facilities in monoclonal antibody (mAb) manufacturing. To help meet these needs, we developed ready-to-use Fibro PrismA units (Fig 1) for capturing mAbs and F_c-containing recombinant proteins. Fibro PrismA units have a protein A cellulose fiber matrix with an open pore structure where mass transfer is governed by convective flow. This structure allows high mAb binding capacities at very short residence times, which results in cycle times of minutes instead of the hours needed for resin-based chromatography. The Fibro PrismA units can be used for up to 200 cycles before disposal, depending on the application.

With HiTrap™ and HiScreen™ Fibro PrismA units connected to ÄKTA™ system you can:

- Purify proteins quickly using rapid cycling chromatography. Cycle times are less than five minutes compared with hours for chromatography resins.
- Achieve high-throughput purification of up to 500 mAbs/wk for clone selection and lead candidate optimization. Each run/purification cycle has real-time UV, pH, and conductivity detection that generates a chromatogram.
- Increase throughput up to 20-fold over that for resin-based chromatography. This cuts weeks from process development lead times.
- Perform a full lifetime study in less than 24 hours.



Fig 1. HiTrap and HiScreen Fibro PrismA units are delivered ready for use.

Product overview

Fibro chromatography, based on electrospun cellulose, offers a large surface area for high binding capacity. The matrix has an open structure with high mechanical strength, which allows high flow rates. Residence times are measured in seconds rather than the minutes required for resin-based chromatography.

The proprietary structure of the cellulose fibers allows the technology to overcome the diffusional and flow limitations of packed bed chromatography purification, as well as the capacity issues of membrane adsorbers and monoliths (Fig 2 and Fig 3). Fibro PrismA has the same engineered PrismA protein A ligand as the MabSelect™ PrismA chromatography resin. This assures optimized binding capacity and excellent alkaline resistance for efficient cleaning-in-place (CIP).







Fig 2. Flow rates in resin-based chromatography (left) are limited by diffusion. The open structure in Fibro fibers (right) allows convective flow and direct mass transfer of the target protein to the ligand immobilized on the fiber surface.



Fibro matrix
Convective mass transfer
Surface area ~ 10 m²/g



Chromatography resins
Binding via diffusion into the particle
Surface area ~ 40 m²/g



Membrane adsorbers
Convective mass transfer
Surface area ~ 0.9 m²/g

Fig 3. Surface area and mass transfer mechanisms for different chromatography base matrices.

HiTrap and HiScreen Fibro PrismA units

These ready-to-use units are designed for research and early process development and are suitable for screening and optimization of process conditions. They are made of medical-grade polypropylene plastic with stoppers at the inlet and outlet.

You can operate both units with either a peristaltic pump or an ÄKTA chromatography system. The HiTrap Fibro unit can also be operated with a syringe and using bidirectional flow. For the smallest possible elution pool volume with the HiScreen Fibro unit, use the defined inlet and outlet to elute proteins in the defined flow direction.

The HiTrap Fibro PrismA unit has one bed, which allows for operation at very high flow rates of 40 matrix volumes (MV)/min on an ÄKTA system. The HiScreen Fibro PrismA unit has two parallel beds that are each twice the thickness of the bed in the HiTrap format (Fig 4). Thus, the recommended flow rate for HiScreen Fibro is lower (8 MV/min) than the flow rate for the HiTrap format.




Fig 4. Cross-sectional image of HiTrap (left) and HiScreen (right) Fibro PrismA units. The bed design in the forthcoming large-scale Fibro PrismA units will be similar to HiScreen Fibro PrismA, but the flow path design will differ.

The characteristics of HiTrap and HiScreen Fibro Prisma units are summarized in Table 1.

Table 1. Main characteristics for HiTrap and HiScreen Fibro Prisma units

| | HiTrap Fibro Prisma | HiScreen Fibro Prisma |
|---|--|----------------------------|
| Matrix | Derivatized electrospun cellulose fibers | |
| Ligand | PrismaA ligand (alkali-stabilized protein A derived from <i>E. coli</i>) | |
| Ligand coupling | Single point attachment | |
| Dynamic binding capacity (DBC) ¹ | ~ 30 mg IgG/mL matrix | |
| Typical DBC/unit ¹ | ~ 12 mg IgG/HiTrap unit | ~ 112 mg IgG/HiScreen unit |
| Cycle time ² | ~ 3 min | ~ 5 min |
| Flow, recommended operating ³ | ≤ 16 mL/min (40 MV/min) | ≤ 30 mL/min (8 MV/min) |
| Matrix volume | 0.4 mL | 3.75 mL |
| Elution volumes ⁴ | ≤ 7 MV | ≤ 4 MV |
| Maximum operating pressure | 1 MPa (10 bar) | |
| Chemical stability | Compatible with aqueous buffers commonly used for protein A chromatography | |
| pH stability, operational ⁵ | 3 to 12 | |
| pH stability, CIP ⁶ | 2 to 14 | |
| Temperature stability, operational | 4°C to 35°C | |
| Temperature stability, storage | 2°C to 8°C | |
| Storage | 20% v/v ethanol in water | |

¹ Determined at 10% breakthrough by frontal analysis in Tris buffer, pH 7.5.

² Purification cycle including steps of: equilibration, sample loading, washing, elution, strip, CIP and reequilibration of the Fibro unit.

³ At room temperature using a buffer with the same viscosity as water.

⁴ HiTrap Fibro and HiScreen Fibro unit housing designs have not been optimized for small elution volumes. But the forthcoming process-scale good manufacturing practices (GMP) compatible Fibro units will be. Elution volumes are expected to be < 3 MV.

⁵ pH range where resin can be operated without significant change in function.

⁶ pH range where resin can be subjected to cleaning- or sanitization-in-place without significant change in function.

The PrismaA ligand

The electrospun cellulose fibers in Fibro are coupled with PrismaA protein A ligand, which is produced in *E. coli*. Fermentation and subsequent purification are performed in the absence of animal products. The ligand has been specifically engineered for enhanced alkali and protease stability. The specificity of binding to the F_c region of IgG is similar to that of the conventional protein A ligand and provides excellent purification in one step. The PrismaA ligand also has affinity for the VH3 chain, so it can be used to purify certain types of antibody fragments.

High binding capacity at very short residence times

The macroporosity and large surface area of the Fibro matrix allow very fast purification with residence times in seconds, not the minutes required with resin-based chromatography (Fig 5). This means that mAbs can be purified up to 20 times faster than with resin-based chromatography. A full mAb purification cycle with equilibration, loading, washes, elution, CIP, and re-equilibration can be performed in minutes rather than hours, as shown in Figure 6.




Fig 5. The immediate mass transfer of the Fibro matrix enables high binding capacity at very short residence times.






Fig 6. Typical purification cycle for high capacity loads using (A) a HiTrap Fibro Prisma unit and (B) a HiTrap 1 mL column containing MabSelect PrismaA resin.

High alkaline stability


With the enhanced alkaline stability of the Prisma protein A ligand, the recommended CIP is 0.5 to 1.0 M NaOH for 30 s to 1 min in every cycle. This means that HiTrap and HiScreen Fibro Prisma units can be confidently cleaned for reuse. If needed Fibro Prisma units can be cleaned with up to 2 M NaOH over hundreds of cycles while still retaining a high binding capacity (Fig 8).




Similar purification performance to resin-based columns



| Product | Recovery (%) |
|-------------------------|--------------|
| HiTrap Fibro PrismA | ~98 |
| HiTrap MabSelect PrismA | ~92 |



| Product | HCP (ppm) |
|-------------------------|-----------|
| HiTrap Fibro PrismA | ~690 |
| HiTrap MabSelect PrismA | ~640 |



A bar chart comparing the leached protein A levels for two chromatography media. The y-axis is labeled 'Leached Protein A (ppm)' and ranges from 0 to 60 in increments of 10. The x-axis has two categories: 'HiTrap Fibro PrismA' and 'HiTrap MabSelect PrismA'. The bar for 'HiTrap Fibro PrismA' is teal and reaches approximately 9 ppm. The bar for 'HiTrap MabSelect PrismA' is dark blue and reaches approximately 24 ppm.

| Media | Leached Protein A (ppm) |
|-------------------------|-------------------------|
| HiTrap Fibro PrismA | ~9 |
| HiTrap MabSelect PrismA | ~24 |

4




Fig 12. Aggregate levels (%) after purification of mAb 1.




Fig 13. Overlaid chromatograms of cycles 1 and 13 of Fibro PrismaA large-scale (0.6 L) prototype unit.

Fibro PrismaA: a scalable protein A fiber chromatography platform

The Fibro technology platform is designed for scalability up to manufacturing scale and compatibility with existing ÄKTA purification systems. The units will range from research scale with a capacity of approximately 12 mg mAb per purification cycle (HiTrap units) to good manufacturing practices (GMP) compatible units that can process up to 2000 L bioreactor harvest, or more than 10 kg mAb in less than 24 hours by cycling multiple times.

Short cycle times enable high-throughput purification of mAbs and opportunities for substantial time savings during research and process development. At larger scales in the biomanufacturing process, operation in a rapid cycling manner allows full utilization of the protein A lifetime during one batch, as well as cost-effective single-use chromatography.

In an independent evaluation, a process-scale Fibro PrismaA prototype unit (0.6 L) was run for 17 cycles on an ÄKTA ready chromatography system equipped with a High Flow kit. The results showed good performance with DBC of 30.6 g/L and cycle times of 7.3 min with 2.8 bar max delta column pressure (dCP) and 4.8 min with 4.2 bar max dCP. Eluate volumes were below 3 MV at > 95% recovery. Buffer usage was 0.65 L/g (~ 18 MV per run). As shown in Figure 13, purification performance was maintained across cycles.

Fibro units are compatible with existing ÄKTA chromatography systems. We recommend using HiTrap Fibro PrismaA with ÄKTA pure/avant 25 and 150. For HiScreen Fibro PrismaA we recommend ÄKTA pure/avant 150. Instrument configurations with predefined UNICORN™ chromatography methods for Fibro are available for these systems. HiTrap and HiScreen Fibro units can also be run on other ÄKTA systems, as shown in Table 2. The forthcoming large-scale GMP compatible Fibro PrismaA units will be compatible with ÄKTA pilot 600 system, ÄKTApurify™ system, and ÄKTA ready system.

Storage

Store HiTrap and HiScreen Fibro PrismaA in 20% ethanol at 2°C to 8°C.

Before use, equilibrate with binding buffer and perform a blank run, including CIP.

Table 2. Fibro units and ÄKTA system compatibility based on recommended flow for each Fibro unit on each system. Blue = preferred system, gray = compatible but not optimal

| System | Chromatography system | | HiTrap Fibro PrismaA (0.4 mL) | HiScreen Fibro PrismaA (3.75 mL) |
|---------------------|-----------------------|-------------------|-------------------------------|----------------------------------|
| | Min flow (mL/min) | Max flow (mL/min) | Recommended flow, mL/min | Recommended flow, mL/min |
| ÄKTA go | 0.01 | 25 | 16 | 25 |
| ÄKTA pure/avant 25 | 0.001 | 25 | 16 | 25 |
| ÄKTA pure/avant 150 | 0.01 | 150 | 16 | 30 |

Ordering information

| Product | Column size | Product code |
|----------------------------|-------------|---|
| HiTrap Fibro Prisma 1 pack | 0.4 mL | 17549855 |
| HiTrap Fibro Prisma 4 pack | 4 × 0.4 mL | 17549856 |
| HiScreen Fibro Prisma | 3.75 mL | Available for orders during Q2. Contact your local sales representative |

[cytiva.com/fibro](https://www.cytiva.com/fibro)

For local office contact information, visit [cytiva.com/contact](https://www.cytiva.com/contact)

Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate. ÄKTA, ÄKTAprocess, HiScreen, HiTrap, MabSelect, and UNICORN are trademarks of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

All goods and services are sold subject to the terms and conditions of sale of the supplying company operating within the Cytiva business. A copy of those terms and conditions is available on request. Contact your local Cytiva representative for the most current information.

CY7145-08Apr20-DF

