
Image Quant TL

(ver.10.2)

Gel and Blot Analysis 簡易操作手順

Windows 版

2023.12

71-4013-32

Image Quant TL 操作手順書

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
The screenshot shows the ImageQuantTL software interface. On the left is a navigation panel with icons for Image, Lanes, BG, Bands, MW, Compare, Normalize, Quant. Cal., Purity, and Results. The main area is divided into three sections: a top-left panel for navigation steps, a large central image display working screen showing a 2D gel image and a 3D surface plot, and a right-side analysis settings panel. The settings panel includes an 'Image Setup' section with instructions on inversion and contrast, an 'Open' button, file information for '1d_Basic.tif', and a highlighted 'Invert all channels' button.

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Large Image display working screen

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The screenshot shows the software interface with the 'Zoom' and 'Fit Image' buttons highlighted. An orange arrow points to the 'Zoom' button with the text '拡大する' (Zoom in). Another orange arrow points to the 'Fit Image' button with the text '元のサイズに戻す' (Return to original size). The main image display shows a 2D gel image with a magnified view of a specific region.

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


The screenshot shows the 'Define Background' dialog box in the software. The 'Background Method' is set to 'Rolling Ball' (①). The 'Radius (Fraction of Lane Length)' is set to 6 (②). The 'Apply' button is highlighted (③). The 'Display Background' checkbox is checked (④). The background shows a gel image with 10 lanes and a corresponding intensity profile graph.

Name	Length	Volume
Lane 1	100.03 mm	413.28133533
Lane 2	100.03 mm	251.84298640
Lane 3	100.03 mm	1032.18839251
Lane 4	100.03 mm	626.17164545
Lane 5	100.03 mm	1181.30407767
Lane 6	100.03 mm	232.70536873
Lane 7	100.03 mm	418.63891987
Lane 8	100.03 mm	182.63100000
Lane 9	100.03 mm	959.15250000
Lane 10	100.03 mm	528.98831216

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

The screenshot shows the 'Select' dialog box with a list of lanes (Lane 1 to Lane 6) and checkboxes. The 'Select' button is highlighted with a red circle and a checkmark, indicating that all lanes are selected. The text '全てのレーンを選択' (Select all lanes) is visible next to the dialog box.

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
The screenshot shows the 'Define Bands' panel with the following settings:

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Normalization Factor = $\frac{\text{Volume of largest overall HK band in reference channel}}{\text{Volume of HK band in each lane in reference channel}}$ > or = 1.0

Normalized Volume = Normalization Factor x Volume of each lane in active channel

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
Normalization Factor = $\frac{\text{Volume of most abundant lane in reference channel}}{\text{Volume of each lane in reference channel}}$ > or = 1.0

Normalized Volume = Normalization Factor x Volume of each lane in active channel

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



The screenshot displays the Cytiva software interface for image analysis. It includes a gel image on the left, a volume profile graph in the center, and a 'Quantity Calibration' window on the right. The calibration window shows a table of known values and a graph of arbitrary units vs volume (%). The main interface also shows a table of detected bands with their positions and volumes.

Lane #	Band #	Position	Dist.	Unit	Type	Volume	Volume - Cali...
L1	1	10.58	mm	detected	60511	5905.6967	
L1	2	22.01	mm	detected	64943	6340.1596	
L1	3	46.40	mm	detected	43226	4211.2087	
L1	4	57.23	mm	detected	16384	1579.8768	
L1	5	69.09	mm	manual	450	16.7029	
L1	6	86.02	mm	detected	40611	3954.9044	
L2	1	11.18	mm	detected	42196	4110.2776	
L2	2	22.69	mm	detected	39651	3860.7125	
L2	3	47.75	mm	detected	27696	2688.8196	
L2	4	57.57	mm	detected	9172	872.9155	
L2	5	86.02	mm	detected	37453	3645.3032	

Quantity Calibration

- Provide known values for individual bands
- Plot known values against measured volumes and fit a curve
- Use the curve to calculate calibrated volumes for all bands in the project

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

The screenshot shows the software interface with the 'Link to Chromatogram' panel on the right. The 'Import Chromatogram' button is highlighted with an orange box. The panel includes a list of steps: 1. Import Chromatogram Curves, 2. Input Sample Information, 3. Add Method Information, 4. Link Chromatogram Peaks to Lanes, 5. Export PDF. Below the steps is a 'Comment' field and a 'Display' section with checkboxes for Image, Lanes, Lanes and Bands, None (image only), Lane Table, Chromatogram, and Lane Profile.

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
The screenshot shows the software interface with the 'Link to Chromatogram' panel on the right. The 'Import Chromatogram' button is highlighted with an orange box. The 'Comment' field is highlighted with an orange box and contains the text 'test'. The 'Lane Table' is highlighted with an orange box and contains the following data:

Lane #	Band #	Volume	Band %	Lane %
L1	1	74.27807143	17.37	16.93
L1	2	84.47850000	19.76	19.26
L1	3	64.03100000	14.98	14.60
L1	4	26.28200000	6.15	5.99
L1	5	178.42900000	41.74	40.69
L2	1	67.47343243	26.98	25.35
L2	2	61.54126316	24.53	23.13
L2	3	51.79548276	20.65	19.46
L2	4	11.79617241	4.70	4.43

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

The screenshot shows the software interface with the following components:

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The screenshot shows the 'Export CSV' dialog box with the following options:

- Detailed Appendices
- Profile
- Profile - Raw
- Volume
- Volume - Raw
- Background
- Export CSV** (button)

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
Field	Value
Image Width	1200
Image Height	1648
Bits per Sample	16
Samples per Pixel	1
Calibrated	no
Calibration Units	n/a
Minimum Value	0
Maximum Value	65535
Photometric Interpretation	0
Comment	Comment = LED&Filter = Chemiluminescence Light = None Filter = Through Tray position = Upper NP lens for tray = No Iris = F0.85 Exposure time = 25.1 s Date time = 2021:03:02 11:17:13 Time zone = Asia/Tokyo(UTC+09:00) Username = Autologin Binning = 5x5 Focus pulse = -3068 Focus mm = 0.0 Software = Amersham ImageQuant 800 Control Software V1.2.0 Firmware1 = 10000 Firmware2 = 01-01 FPGA = 8 Serial number = 16120116 Instrument = Amersham ImageQuant 800 Platform = Flat_chem_upperFlatframe date = 2019:08:28 12:42:04 Invert = No A/RGB = No OD parameter = OD Conversion coefficient = 1/1.4 OD mode = NA Orientation = top-left Image split = none Range low = 64 Range high = 25536 Imaging mode = Chemiluminescence Method = Exposure mode = Auto Image Area = 80x110 Exposure Number = No Saturated in JPEG = Highlight Dynamicrange = ON Setting1 = 1251 Setting2 = No Setting3 = 672 Setting4 = 2251 Correction1_1 = 2021:01:07 11:38:34 Correction1_2 = 2021:01:07 11:38:34 Correction2_1 = 18411477 Correction2_2 = 13774826 Correction2_3 = 143114887 Correction3 = 2020:01:27 16:05:54 Correction4 = 963 Correction5_1 = No Correction5_2 = No Correction6 = 2.5
Imager Make	Unknown
Acquisition Date	Unknown
Imager Model	Amersham ImageQuant 800
Acquisition Time	Unknown
Resolution	143.4720 dots / cm
Pixel Size	69.7 µm

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Trademarks

For local office contact information, visit cytiva.com/contact

Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate.

Amersham and ImageQuant are trademarks of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

Apple is a trademark is a trademark of Apple Incorporated. Windows is a trademark of Microsoft Corporation.

All other third-party trademarks are the property of their respective owners.

© 2023 Cytiva

All goods and services are sold subject to the terms and conditions of sale of the supplying company operating within the Cytiva business. A copy of those terms and conditions is available on request.

Contact your local Cytiva representative for the most current information.

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