
Cytiva Webinar

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Biacore SPR system

Single needle to eight needle parallel SPR analysis instruments

Structure –function analysis, binding, affinity

Biacore software

System-specific software and Biacore™ Insight Evaluation Software with software extensions

Screening, titer, affinity, kinetics, Sensorgram Comparison

Biacore consumables and sensor chips

Save time and effort with functionalized surfaces, ready-to use capture kits and sensor surfaces for a broad range of molecules

Characterization and quality control, epitope binning, kinetics, concentration, Fc-receptor binding

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IQ/OQ/IPQ = installation qualification/
operational qualification/initial performance qualification

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Assessment of protein concentration using a calibration-free method

Kevin Lindqvist, Åsa Sparro, Robert Karlsson, Anna Sylvan, Anita Larsson, and Bwa Pol. GE Healthcare, Uppsala, Sweden.

Introduction

In development and production of biopharmaceuticals it is important to determine the concentration of functional glycoprotein molecules in a given sample. Many available methods for protein concentration measurement do not distinguish between active and inactive molecules, or rely on comparison with a standard that may or may not be a mimic of the actual sample. The method discussed here, Calibration-Free Concentration Analysis (CFCA), provides a highly sensitive and robust bio-sensor-based determination of the protein concentration and does not require a standard. The method relies on measurement of analyte binding to a target immobilized on a sensor surface at varying flow rates, under conditions where the observed rate of binding is partially or completely limited by transport of analyte molecules to the sensor surface. This transport is a diffusion-controlled process that depends on the analyte concentration. The concentration is obtained by running Biacore® binding experiments at different flow rates and fitting the binding data to the model describing the process.

In factor in instrument samples are injected in a micro flow system and transported in laminar flow to the sensor surface. Molecules reach the sensor surface from bulk solution by a diffusion-controlled transport process. In addition to the concentration of analyte molecules, factors influencing the transport rate include the diffusion coefficient, flow cell dimensions, and flow rate. The balance between their transport rate and the binding rate determines whether the observed binding will be transport limited or reaction limited. For successful CFCA, the observed binding rate must be at least partially limited by transport. In this paper, we discuss how suitable flow rates, immobilization levels, and binding data are selected for successful implementation of CFCA.

Immobilization level and concentration range

Calibration-free concentration analysis of 2-fold dilution series of β_2 -microglobulin was performed on four surfaces immobilized to increasing levels of anti- β_2 -microglobulin in a rebinding. Immobilization levels are expressed as response units (RU) divided by the molecular weight of interaction partner, M_w, to compensate for the difference in response between differently sized molecules. The protein solutions were injected at two flow rates, 5 and 90 μ l/min.

The concentration determined by the method is shown in each panel. Sensorgrams from different flow rates are obtained under equivalent conditions that return the same calculated concentration after adjustment for the dilution factor. The concentration given by sensorgrams from red and green sensorgrams are compared.

The most apparent visual indicator of data suitable for CFCA is a significant separation of the binding curves obtained at high and low flow rates.

Flow rates

Concentrations were determined for two independent antibody model systems using binding data obtained at flow rates of 5, 10, 20, 40, 60, 80 and 100 μ l/min. The binding data was analyzed using 23 pairwise combinations of flow rates, as well as in the whole data set of 7 flow rates. Results were normalized to the value obtained from 7 flow rates (shown as red line with error lines to indicate ± 1 , SD) to allow comparison between the model systems.

The results indicate that reliable concentration data is obtained when using two flow rates that differ by a factor of 10 or more.

Concentration of 8 commercial antibodies

The CFCA of eight commercial antibodies, denoted as A to H, was performed under transport limited conditions using anti-Fc antibody and protein A. The concentration values were compared to the concentrations determined by manufacturers. Activity is given as the ratio of measured and specified concentration, respectively. The concentrations of all antibodies were in good agreement when measured against two different interaction partners (although could differ as much as 50% when compared to the concentration specified by manufacturers).

How to recognize transport limited data?



Sensorgrams obtained at two flow rates.

- Transport limited data:
 - High binding at low flow rate
 - Low binding at high flow rate
 - Highly sensitive to flow rate
 - Not suitable for CFCA
- Reaction limited data:
 - High binding at high flow rate
 - Low binding at low flow rate
 - Not sensitive to flow rate
 - Not suitable for CFCA
- Transport and reaction limited data:
 - High binding at low flow rate
 - Low binding at high flow rate
 - Highly sensitive to flow rate
 - Highly sensitive to flow rate
 - Not suitable for CFCA

Conclusions

- Calibration-free concentration analysis on Biacore systems opens new possibilities for concentration measurements that relate directly to determination of the concentration of a particular binding site.
- Measurements are direct, label-free, and do not require a calibration standard.
- Robust and reliable results are obtained when experiments are performed using:
 - at least two flow rates, where the ratio between flow rates is > 10
 - a high surface concentration, typically > 25 RU/M_w (immobilized level/M_w of interaction partner)
 - protein dilutions that give concentrations \leq 30 nM
- In the example, the concentrations of eight antibodies were in good agreement when measured against two different interaction partners targeting the antibody Fc part. However, concentrations determined by CFCA differed by up to 50% from the values specified, depending on the method used. More work is required to explain the discrepancies.

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



High-resolution SPR-based surrogate potency assay to facilitate comparability and biosimilar studies
 Asa Frostell, Veronica Fridh, and Robert Karlsson
 GE Healthcare Bio-Sciences AB, Arturqvist 70, 751 84 Örebro, Sweden

Introduction
 An assay platform is needed to measure and monitor critical quality attributes (CQAs) for antibody development. The present study describes a high-resolution SPR-based surrogate potency assay for full-scale comparability and stability studies.

Surrogate potency assay using a biotinylation-based protein-coupled antibody
 Biotinylation of antibodies was used to facilitate the immobilization of antibodies on a sensor chip. The resulting biotinylated antibodies were used to immobilize antibodies on a sensor chip.


Conclusions
 A high-resolution SPR-based surrogate potency assay was developed. The assay was used to measure and monitor critical quality attributes (CQAs) for antibody development. The assay was used to measure and monitor critical quality attributes (CQAs) for antibody development.

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


1) $K_D/k_d/k_a$
1:1 binding data


2) Report points
A limited subset of data

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

Immobilization

Regeneration, if needed


Response data from samples and reference substance

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Immobilization levels of different lots of Protein A from one supplier

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


Figure 2: Optimization of ligand immobilization conditions. Maximal binding of MT201 to protein A was achieved using 10 µg/ml protein A at a flow rate of 5 µl/minute for 7 minutes.

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

Change of standard curves measured on Sensor Chip CM5 over 2 months storage in HBS-EP buffer

Application Note 38

Biacore systems

Rapid development of a GMP-compliant assay for the determination of antibody concentration

- Easy, rapid validation of assay performance
- Assays compliant with current ICH guidelines

Introduction

Manufacturing regulations require that the concentration of proteins produced by recombinant or hybridoma technology is monitored throughout the production process using validated assays. Here, we describe the process leading to the establishment of a validated assay for determining the concentration of a monoclonal, fully human anti-*EpCAM* IgG1 antibody (Naundorf et al., 2002) intended for the treatment of prostate cancer, from harvest samples, process controls and final product. The conditions were optimized within three months to deliver an accurate and reliable assay.

Assay validation: preliminary conditions

Although no definitive guidelines on the validation of binding data are available today, the FDA documents, *Guidance for Industry: Bioanalytical Method Validation and Guidance for Industry: Analytical Procedures and Methods Validation Chemistry, Manufacturing, and Controls Documentation* provide generally accepted guidelines for binding assay design. Until consensus is established, individual QA departments are required to comply with the current ICH guidelines for assay validation.

Once the basic design of the assay is decided, the precise conditions under which it will be run must be determined empirically. Here, an assay was constructed according to the FDA guidelines above, in which a humanized monoclonal antibody is captured on Sensor Chip CM5 after initial immobilization of protein A serving as the IgG-specific binding molecule. Immobilization, surface regeneration, and sensor

Immobilization levels of different lots of protein A from one supplier


Figure 1. Immobilization levels of 16 lots of protein A from a single supplier on certified Sensor Chip CM5. The mean immobilization level using 10 µg/ml protein A at a flow rate of 5 µl/minute for 7 minutes was 1970 ± 137 RU (CV 7.0%).

chip storage were all investigated and the optimal conditions for the final assay were determined from these data. All experiments were performed using certified Sensor Chip CM5 and a Biacore™ C system.

Optimization of protein A immobilization conditions

Protein A was immobilized on certified Sensor Chip CM5 in order to measure the consistency of immobilization levels between different batches. Sixteen lots of protein A prepared at a concentration of 10 µg/ml and injected over the sensor chip at a flow rate of 5 µl/minute for 7 minutes gave immobilization levels of between 1500 and 2500 RU (Figure 1). There were no significant differences in the immobilization levels of different lots of protein A under these conditions. It is recommended that routines are established to achieve a similar level of protein A immobilization before further assay optimization.

Acknowledgement

We gratefully acknowledge Birgit Vetter, Susanne Schade, Dr Thomas Urbig and Dr Frank Handlham from Micromet AG, Munich for providing the experimental data for this Application Note.

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

Figure 3. Test of receptor stability by repeated injection of analyte.

Application Note 48 Biacore systems

Validation of a concentration assay using Biacore C

- Guideline for development of a GxP-compliant concentration assay
- Support for informed decision making in drug development

Introduction

In the process of measuring active concentrations of pharmaceuticals, high demands are placed on the quality of the data on which any decision to proceed to clinical trials will be based. Data gathered using any technology must be accurate and precise. Method validation (defining the acceptable limits within which data is deemed reliable) is therefore a crucial part of the overall process through which a prospective drug passes before it is tested in clinical trials. Samples for pharmaceutical concentration measurements are also taken from reactor runs, to monitor the development of purification processes, in stability studies and in hybridoma screening. Here we present the results obtained during a validation process within the pharmaceutical industry for a concentration assay using Biacore™ C system. The results confirm that the quality of performance of this assay was sufficiently high to meet GxP demands placed on drug companies by regulatory authorities.

Biacore C is specifically designed for concentration analysis of pharmaceuticals. The instrument software supports GxP-regulated working processes with user-friendly wizards for assay development, immobilization strategies and data analysis. In addition, as current regulations demand restricted access to data generated by analytical instruments, Biacore C software makes it possible to lock and fix assay templates.

In the example cited here a receptor/ligand system was investigated. Both receptor and ligand are proteins. The receptor protein was immobilized at high density on the surface of an activated sensor chip. Reference material,

quality control (QC) samples and test samples were then diluted and injected over the receptor surface to generate sensorgrams in which binding data is recorded over real time. The initial binding rate of the ligand to immobilized receptor (sensorgram slope) was measured under conditions of limiting mass transfer (diffusion). The rate, expressed as resonance units per second (RU/s) correlated directly with active ligand concentration. The binding rates of the samples were converted into active concentration using the reference standard curve where the binding rate of reference material was plotted against the concentration. The final results were expressed either as percent binding of sample relative to reference material or as a concentration.

Concentration assay validation: the cornerstones

Several criteria were tested over the course of the validation outlined below. Experiments should include, but are not necessarily limited to the following tests:

Parameter	Question addressed
Specificity	To what extent does the assay cross-react with related targets?
Linearity, range, detection and quantitation limits	To what degree of certainty can concentration measurements be integrated on a standard curve and over what range can concentration measurements be confidently stated?
Accuracy	Do the determined concentration measurements lie within acceptable error limits when compared with nominal concentrations?
Precision	Can similar concentrations be measured with an acceptable level of reproducibility between users on different occasions (between assay variations and within the same assay/assay variation)?
Robustness	Is the assay robust enough to perform under sub-optimal conditions, e.g. after repeated freeze/thaw cycles of the sample or if an alternative buffer system is used?


Acknowledgement

We gratefully acknowledge Dr Ralph Abraham and Mary Quimby of Bristol-Myers Squibb Inc., Pharmaceutical Research Institute, Bi analytical Sciences, Hopewell, New Jersey, USA for providing the data for this Application Note.

Biacore



Affinity / Kinetics assay development

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

Rate constants, level of robustness:

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Kinetics

Parameter	Acceptance criteria, an example
k_a ($M^{-1} s^{-1}$)	5×10^5 to 1×10^6
k_d (s^{-1})	2×10^{-3} to 4×10^{-3}
R_{max} (RU)	60 to 80

Affinity



Parameter	Acceptance criteria, an example
K_D (M)	1×10^6 to 2×10^6
R_{max} (RU)	50 to 100

Additional criteria

Parameter	Qualification Target
Immobilization level for further ligand capture	% CV \leq 40%
Immobilization/capture levels for analyte binding	% CV \leq 30%
Chi ² for kinetic and affinity fit	The square root of Chi ² for the fit should be insignificant, < 5% to 10%, in relation to the binding level

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5

Case study

Validation of kinetic assays using Biacore X100

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

Acceptance criteria for kinetic based on statistical results from the validation

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
Insulin binding to anti-insulin R-INS-6

- k_a (1/Ms): $> 8.56E+06$
- k_d (1/s): < 0.118


	k_a (1/Ms)	k_d (1/s)	K_D (M)
Run 1	1.55E+07	0.09408	6.09E-09
Run 2	1.19E+07	0.07399	6.22E-09
Run 3	1.75E+07	0.09741	5.58E-09
Run 4	1.61E+07	0.09405	5.86E-09
Run 5*	(8.01E+06)	(0.1281)	(1.60E-08)
Run 6	1.38E+07	0.07833	5.68E-09
Run 7	1.25E+07	0.06847	5.49E-09
Mean	1.45E+07	8.44E-02	5.82E-09
CV (%)	13.67	13.29	4.55

*Not approved due to high U-value

Insulin binding to anti-insulin R-INS-5

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Anti-insulin R-INS-6	k_a (1/Ms)	k_d (1/s)
Acceptance criteria	$> 8.56E+06$	< 0.0118
New lot		
Measured value	$1.94E+07$	0.01094

Anti-insulin R-INS-5	k_a (1/Ms)	k_d (1/s)
Acceptance criteria	$> 6.84E+05$	$< 7.38E-03$
New lot		
Measured value	$1.14E+06$	$6.64E-03$

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Curve comparison for QC

Karen Duus, Research Scientist, Novo Nordisk

Using SPR for FDA Approved Batch Release and Assay Validation.

Matthew Rutter, QC Analyst, GSK

Global Analytics Update: Regulatory, QC and Validation White Paper Review

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