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Medium	ActiCHO P with 6 mM L-glutamine
Culture chamber	50 L Cellbag bioreactor
Rocking	
Low-volume stage	20 rpm, 6° rocking angle
Full-volume stage	20 rpm, 8° rocking angle
Gas flow	0.3 L/min
Seed cell concentration	0.3×10^6 viable cells/mL
pH set point	
Low-volume stage	No pH control (initial 7.5% CO ₂ in headspace)
Full-volume stage	pH controlled at 7.0 with CO ₂
Dissolved oxygen (DO) set point	40%
Working volume	2.6 L for initial 3 d in shaker flasks, 15 L to target cell concentration in ReadyToProcess WAVE 25
Target cell concentration at transfer	3.5 to 4.5×10^6 viable cells/mL
Viability criteria at transfer	> 95%

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	XDR-200	WAVE 25
Medium	ActiCHO P with 6 mM L-glutamine	ActiCHO P with 6 mM L-glutamine
Culture chamber	Standard bag assembly: Development	20 L Cellbag bioreactor with DO and pH sensors
Starting viable cell concentration	0.35 to 0.45×10^6 viable cells/mL	0.35 to 0.45×10^6 viable cells/mL
Inoculum volume	15 L	7.0 L transferred from XDR-200
Medium fill volume	130 L	N/A
Starting operating volume	145 L 138 L after transfer to satellite bioreactor	7.0 L
Impeller/rocker speed	150 rpm	20 to 30 rpm, 8° rocking angle, 30% rocking motion
pH set point	7.1	7.1
pH PID* settings	P: 0.5 I: 1.0 D: 0	Factory settings (automatically computed by the system)
Base for pH control	7.5% w/v NaHCO ₃	7.5% w/v NaHCO ₃
DO set point	40%	40%
DO PID* settings	P: 0.254 I: 4.083 D: 0	Factory settings (automatically computed by the system)
Air gassing	2 L/min head-space gassing DO-controlled gassing administered through 1 mm sparger Controller output: 0% to 20%: 1 to 5 L/min (sparged) 20% to 100%: 5 L/min (sparged)	On demand (from DO and pH control) 0.2 L/min total gas flow
Oxygen flow	On demand (from DO PID control) administered through 20 μm sparger Controller output: 0% to 20%: 0 L/min 20% to 100%: 0 to 5 L/min (sparged)	On demand (from DO PID control) max. 50% 0.2 L/min total gas flow
Antifoam	Medical antifoam C (Sigma-Aldrich) (added as required)	Medical antifoam C (Sigma-Aldrich) (added as required)
Feeds	ActiCHO Feed-A ActiCHO Feed-B 400 g/L D-glucose	ActiCHO Feed-A ActiCHO Feed-B 400 g/L D-glucose
Feed starting point	Day 2	Day 2
Daily feed volumes	ActiCHO Feed-A: 5.915 L ActiCHO Feed-B: 0.592 L Glucose: as required with target of 6 g/L excluding ActiCHO Feed-A	ActiCHO Feed-A: 0.295 L ActiCHO Feed-B: 0.0295 L Glucose: as required with target of 6 g/L excluding ActiCHO Feed-A
Harvest criteria	13 elapsed days and/or when culture viability was < 60%	13 elapsed days and/or when culture viability was < 60%

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Measured parameter	Assay type	Type of sample
Cell concentration and viability	Bioprofile FLEX analyzer (Nova Biomedical Corp.)	Cell suspension
Glucose	Bioprofile FLEX analyzer	Cell suspension
Lactate	Bioprofile FLEX analyzer	Cell suspension
Ammonium	Bioprofile FLEX analyzer	Cell suspension
Glutamine	Bioprofile FLEX analyzer	Cell suspension
Glutamate	Bioprofile FLEX analyzer	Cell suspension
Osmolality	Bioprofile FLEX analyzer	Cell suspension
Product concentration	CEDEX Bio analyzer (Roche)	Clarified cell culture supernatant

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