
Biacore 1 series

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

3. Insight Control Softwareの起動・ログイン

A screenshot of the Biacore Insight Control Software login interface. The screen displays the Cytiva logo and the text "Biacore™ Insight Control Software Version 5.0.16.21762". There are input fields for "User id" (containing "XXXXXX / 12345"), "Domain" (containing "XXXXXX"), and "Password". Below these fields, there are options for "Database" (set to "Local Biacore Insight database") and "Selected extensions" (listing "Biacore Intelligent Analysis™", "Concentration & Potency", "Data Integration", "Epitope Binning", and "Extended Screening"). At the bottom, there are "Log in" and "Cancel" buttons.

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
1-3. Biacore Insight Control Software

Instrument control画面

The screenshot shows the Biacore Insight Control Software interface. The top navigation bar includes the Cytiva logo, the software name, user information, and a help menu. Below the navigation bar are four main panels, each highlighted with an orange border:

- Activity queue:** Contains a 'Change solutions' task with a 6-minute progress bar and an 'Interactive run' button. An 'Add activity' button is at the bottom.
- System setup:** Lists tasks such as 'Set flow cell temperature', 'Set sample compartment temperature', 'Change chip', 'Change solutions', 'Immobilization checkpoint', and 'Wait'.
- Maintenance:** Lists tasks such as 'Desorb', 'Desorb and sanitize', 'System check', 'Normalize', and 'Shutdown'.
- Start:** Contains 'Interactive run' and 'Method' buttons.

Below the panels, there are four descriptive labels in Japanese:

- Activity queue:** 各種コマンドの予約 (Reservation of various commands)
- System setup:** 温度設定 (Temperature setting), センサーチップのドック (Sensor chip docking), バッファ置換など (Buffer replacement, etc.)
- Maintenance:** Desorb (週一) (Desorb (weekly)), Desorb and Sanitize (月一) (Desorb and Sanitize (monthly)), System Check (適宜) (System Check (as appropriate))
- Instrument states:** システムの状態 (System status)

At the bottom, a status bar displays system information: 'Running Change solutions...', 'Current buffer Buffer', 'CAP: 3/3/2023 6:09:55 PM', 'chip information', 'Flow cell 25.0 °C (25)', 'Sample compartment 25.0 °C (25)', 'Illumination on/off', and 'Hotel door is closed'.

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
The screenshot shows the Biacore Insight Control Software interface. The top navigation bar includes the Cytiva logo, the software name, a user profile, and a help icon. Below this is a secondary navigation bar with tabs for "Instrument control", "Methods", "Runs", and "Action history". The "Methods" tab is active. In the "Methods" section, there are "Open" and "New" buttons, with "New" highlighted by a red box. The main content area is divided into two sections: "Biacore 8 series" and "Biacore 1 series". Under "Biacore 1 series", a tree view shows "Kinetics / affinity" selected. To the right, a table lists various method templates with their names and descriptions. An "Item details" panel is visible on the right side of the table.

Name	Description
Kinetic screen	Single concentration injections for rapid estimation of kinetic properties. Includes...
Single-cycle kinetics	Single-cycle kinetics for protein samples.
Multi-cycle kinetics/affinity	Multi-cycle kinetics/affinity for protein samples. Includes regeneration.
Single-cycle kinetics using capture	Single-cycle kinetics for protein samples using captured ligand.
Multi-cycle kinetics/affinity using capture	Multi-cycle kinetics/affinity for protein samples, using captured ligand. Includes...
Fast kinetics 10 Hz	Multi-cycle kinetics for fast dissociating analytes. Using data collection rate 10
Fast kinetics 40 Hz	Multi-cycle kinetics for fast dissociating analytes. Using data collection rate 40

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Property	Variable	Value	Flow path	Predip
Solution	<input type="checkbox"/>	Capture solution	2	<input checked="" type="checkbox"/>
Contact time	<input type="checkbox"/>	60 s		
Flow rate	<input type="checkbox"/>	10 µl/min		
Concentration	<input type="checkbox"/>	nM		
Molecular weight	<input type="checkbox"/>	Da		

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
Fig 4.1. Ligand is concentrated on the surface through electrostatic attraction when the pH lies between the isoelectric point of the ligand and the pK_a of the surface. If the pH is too low or too high, ligand will not be concentrated on the surface.

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

Instrument control | Methods | Runs | Action history | Preferences

Immobilization | Open | New

Immobilization - 1 Series | 1. Method definition | 2. Positioning and plate layout | Send to queue | Compatibility

Bottles: Buffer bottle 200 ml Buffer, Water bottle 200 ml Water

Sort positions: Ascending, Descending | Estimated run time 32 min

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
Cycle	Channel	Flow cell	Sensorgram type	Analysis step purpose	Analysis step name	Excluded	Curve markers	Capture 1 Solution	Analyte 1 Solution	Analyte stability early_1 Relative (RU)
8	5	2-1	Reference subtracted	Analysis	Analysis 1			Ligand 3	beta-2micro	2.4
8	6	2-1	Reference subtracted	Analysis	Analysis 1			Ligand 3	beta-2micro	2.6
8	7	2-1	Reference subtracted	Analysis	Analysis 1			Ligand 3	beta-2micro	2.6
8	8	2-1	Reference subtracted	Analysis	Analysis 1			Ligand 3	beta-2micro	2.3
9	8	2-1	Reference subtracted	Analysis	Analysis 1			Ligand 3	beta-2micro	30.7
9	7	2-1	Reference subtracted	Analysis	Analysis 1			Ligand 3	beta-2micro	29.4
9	5	2-1	Reference subtracted	Analysis	Analysis 1			Ligand 3	beta-2micro	28.1
9	6	2-1	Reference subtracted	Analysis	Analysis 1			Ligand 3	beta-2micro	28.2

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
4-1. 抗体のKinetics screen測定

Control Softwareの構造～Method definitions

The screenshot shows the 'Method Builder - 1 Series' interface. At the top, there are tabs for '1. Method definition', '2. Variables and positioning', '3. Cycle overview', and '4. Plate layout'. Below these, a 'Startup' and 'Analysis' step are visible. The 'Analysis' step is expanded to show three sub-steps: 'Capture 1', 'Analyte 1', and 'Regeneration 1'. A table at the bottom right explains these steps.

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4-1. 抗体のKinetics screen測定

Variables and positioning

The screenshot displays the 'Method Builder - 1 Series' interface. The '2. Variables and positioning' step is highlighted. The main area shows a table for 'Capture 1 - Solution' and 'Analyte 1 - Concentration (nM)'. The 'Capture 1 - Solution' column has a dropdown menu (labeled 1) and the 'Analyte 1 - Concentration (nM)' column has a dropdown menu (labeled 2). To the right, there are 'Positioning settings' for 'Plate 1' and 'Rack 1', including a grid of wells (labeled 3) and a 'Positioning settings' sidebar (labeled 4).

No	Capture 1 - Solution	Control	Analyte 1 - Concentration (nM)
1	Sample 1	▼	0
2	Sample 1	▼	1
3	Sample 2	▼	0
4	Sample 2	▼	1
5	Sample 3	▼	0
6	Sample 3	▼	1
7	Sample 4	▼	0
8	Sample 4	▼	1
9	Sample 5	▼	0
10	Sample 5	▼	1
11	Sample 6	▼	0
12	Sample 6	▼	1
13	Sample 7	▼	0
14	Sample 7	▼	1
15	Sample 8	▼	0
16	Sample 8	▼	1
17	Sample 9	▼	0

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
2. Send to queueで測定開始

The screenshot shows the 'Kinetic screen' software interface. The top navigation bar includes 'Instrument control', 'Methods', 'Runs', 'Action history', and 'Preferences'. Below this, there are buttons for 'Kinetic screen', 'Open', and 'New'. The 'Method Builder - 1 Series' section has four steps: '1. Method definition', '2. Variables and positioning', '3. Cycle overview', and '4. Plate layout'. The '4. Plate layout' step is highlighted with an orange box and a circled '3'. A 'Send to queue' button is also highlighted with an orange box and a circled '3'. The 'Plate 1' grid shows a 12x8 layout with columns A-H and rows 1-12. A table on the right explains the interface elements with circled numbers 1-3.

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


The screenshot displays the Biacore Insight Evaluation Software interface. The top navigation bar includes the Cytiva logo, the software name, and user information. Below this, a menu bar allows switching between different evaluation modes, with 'Evaluation - Kinetic screen' selected. A 'Thumbnails' section shows a grid of small plots for various ligand conditions. A green arrow points from the selected thumbnail to a larger 'Sensorgram' plot. Below the sensorgram, the 'Fit details' section is visible, showing parameter values for the selected fit.

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Biocore Insight Evaluation Software - Kinetic screen large

cytiva Biacore™ Insight Evaluation Software User HELP

Create new evaluation Open existing evaluation Evaluation - Kinetic screen large Action history

QC - Sensorgram QC - Capture baseline QC - Binding to reference Evaluation - Capture level Evaluation - Kinetic screen Home

Settings

Data grouping (32) Serial Parallel/2D

Injection assignment

Kinetics/Affinity mode Settings apply to selected series

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
A screenshot of the software interface for 'Interactive run' configuration. The interface is divided into several sections: 'Activity queue' on the left, 'Checklist' in the top middle, 'Bottles' in the middle, 'Data collection rate' below that, 'Use flow cells / Reference' below that, 'Concentration units' below that, and 'Tray setup' on the right. A 'Ready to start' button is at the bottom. Numbered callouts (1-6) are placed over the interface to correspond to the table above: 1 points to the 'Use buffer position' dropdown (A, B, C, D); 2 points to the 'Data collection rate' dropdown (1 Hz); 3 points to the 'Use flow cells / Reference' dropdowns (1,2 and 1); 4 points to the 'Concentration units' 'Change units' button; 5 points to the 'Tray setup' grid; 6 points to the 'Ready to start' button.

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The screenshot shows the 'Interactive run' interface for a Cytiva instrument. At the top, there are tabs for 'Instrument control', 'Methods', 'Runs', and 'Action history', along with a 'Preferences' gear icon. Below this, the current run is identified as 'Interactive run 5/10/2023 5:29:04 PM'. There are buttons for 'Sensorgrams' and 'Run properties', and 'End run' and 'Abort run' buttons. A list of flow cells (Flow cell 1 to 6) is shown, with 'Flow cell 1 Ref' and 'Flow cell 2' selected. A '1' is circled next to 'Flow cell 1 Ref'. Below the flow cells is an 'Add command' dropdown menu, circled with a '3'. To the right of the flow cells is an 'Add cycle' dropdown menu, circled with a '4'. On the left, there is a 'Settings' section with a 'Sensorgrams' subsection where 'Flow cell 1' and 'Flow cell 2' are checked, with a '2' circled next to 'Flow cell 2'. Below the settings are buttons for 'Command sequence', 'Sensorgrams', 'Cycle event log', 'Injection pause controls', and 'Response table'. The main area displays a sensorgram graph titled 'Analysis 1' with 'Response' on the y-axis and 'Time' on the x-axis. The graph shows a baseline at 0.0 with a small peak at the start. A legend indicates 'Fc=1' (orange), 'Fc=2' (pink), and 'Fc=2-1' (green).

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
5-2. Kinetics/Affinity (K_D , k_a , k_d)測定

Instrument control Methods Runs Action history Preferences

Single-cycle kinetics using capture Open New

Method Builder - 1 Series

1. Method definition 2. Variables and positioning 3. Cycle overview 4. Plate layout Send to queue

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Variables and positioning

Method Builder - 1 Series

1. Method definition 2. Variables and positioning 3. Cycle overview 4. Plate layout

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

1. Select runs 2. Select evaluation method

Search...

Select runs from multiple folders Show 20 items 1-16

- Root
 - Demo Data
 - Biacore 8K Demo data
 - Biacore8K Training data**
 - IE Concentration Demo Files
 - IE Epitope Binning Demo Files
 - IE Intelligent Analysis Demo Files

Name
System check 6/9/2017 10:29:35 AM
SCK demo 11/30/2016 1:23:47 PM
Fragment affinity screen 11/2/2016 3:33:31 PM
Fragment clean screen 11/2/2016 12:31:04 PM
Immobilization Clean Screen 11/2/2016 11:09:57 AM
Single-cycle kinetics 011116 11/1/2016 3:40:51 PM

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1. Select runs 2. Select evaluation method

User defined **Predefined**

- Empty method
- Surface preparation
- Assay development
- Binding screen
- Kinetics / affinity
 - Fragment
 - LMW
 - Antibody / general**

Name
Kinetic screen - Evaluation method
Single-cycle kinetics - Evaluation method
Multi-cycle kinetics - Evaluation method
Multi-cycle affinity - Evaluation method
Single-cycle kinetics using capture - Evaluation method
Multi-cycle kinetics using capture - Evaluation method
Multi-cycle affinity using capture - Evaluation method

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The screenshot displays the Biacore Insight Evaluation Software interface for a single-cycle kinetics analysis. The main window is titled "Evaluation - Single-cycle kinetics 011116 11/1/2016 3:40:51 PM" and shows a "QC - Binding to reference" view. The interface is divided into several sections:

- Settings (left sidebar):** Includes options for "Data grouping (8)", "Serial/Parallel/2D", "Injection assignment", "Kinetics/Affinity mode" (set to Kinetics), "Blank settings", "Fit models" (Perform fit [0]), "Initial values", and "Quality prediction: None".
- Thumbnails (top right):** A row of seven small plots showing sensorgrams and fits for different analyte/ligand ratios (1:1 and 2:1). One thumbnail is highlighted with a red circle.
- Sensorgrams (center):** A large plot titled "Sensorgrams" showing "Relative response" vs "Time" for "Analyte; Ligand 1; 1:1 binding". The plot shows a typical binding curve with a fit line. Below the plot, parameters are listed: "1:1 binding; ka=2.58e+05; kd=7.99e-04; Rmax=733.2; KD=3.10e-09".
- Fit detail (bottom right):** A table showing "1:1 binding Sensorgram parameters" with columns for "Curve", "Conc (M)", and "f (µl/min)". The table lists concentrations from 2E-09 to 3.2E-08 M and a response of 30 µl/min.

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