
Biacore 1 series

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

3. Insight Control Softwareの起動・ログイン

A screenshot of the software's login interface. At the top, it says "cytiva" and "Biacore™ Insight Control Software Version 5.0.16.21762". There are input fields for "User id" (containing "XXXXXX / 12345"), "Domain" (containing "XXXXXX"), and "Password". Below these are options for "Database" (set to "Local Biacore Insight database") and "Selected extensions" (listing "Biacore Intelligent Analysis™", "Concentration & Potency", "Data Integration", "Epitope Binning", and "Extended Screening"). At the bottom, there are "Log in" and "Cancel" buttons.

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
1-3. Biacore Insight Control Software

Instrument control画面

The screenshot shows the Biacore Insight Control Software interface. The top navigation bar includes the Cytiva logo, the software name, user information, and a help menu. Below the navigation bar are four main panels, each highlighted with an orange border:

- Activity queue:** Contains a 'Change solutions' task with a 6-minute progress bar and an 'Interactive run' task. An 'Add activity' button is at the bottom.
- System setup:** Lists tasks such as 'Set flow cell temperature', 'Set sample compartment temperature', 'Change chip', 'Change solutions', 'Immobilization checkpoint', and 'Wait'.
- Maintenance:** Lists tasks such as 'Desorb', 'Desorb and sanitize', 'System check', 'Normalize', and 'Shutdown'.
- Start:** Contains 'Interactive run' and 'Method' buttons.

Below the panels, there are four descriptive labels in green text:

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At the bottom of the interface, a status bar displays various system metrics and indicators, including 'Running Change solutions...', 'Current buffer Buffer', 'CAP: 3/3/2023 6:09:55 PM', 'chip information', 'Flow cell 25.0 °C (25)', 'Sample compartment 25.0 °C (25)', 'Illumination on/off', and 'Hotel door is closed'.

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
The screenshot shows the Biacore Insight Control Software interface. The 'Methods' tab is selected in the top navigation bar. Below the navigation bar, there are 'Open' and 'New' buttons. The 'New' button is highlighted with a red box. The main content area is divided into two sections: 'Biacore 8 series' and 'Biacore 1 series'. Under 'Biacore 1 series', the 'Kinetics / affinity' folder is expanded, showing several sub-methods. A table lists these methods with their names and descriptions. The 'Item details' panel on the right is empty.

Name	Description
Kinetic screen	Single concentration injections for rapid estimation of kinetic properties. Includes...
Single-cycle kinetics	Single-cycle kinetics for protein samples.
Multi-cycle kinetics/affinity	Multi-cycle kinetics/affinity for protein samples. Includes regeneration.
Single-cycle kinetics using capture	Single-cycle kinetics for protein samples using captured ligand.
Multi-cycle kinetics/affinity using capture	Multi-cycle kinetics/affinity for protein samples, using captured ligand. Includes...
Fast kinetics 10 Hz	Multi-cycle kinetics for fast dissociating analytes. Using data collection rate 10
Fast kinetics 40 Hz	Multi-cycle kinetics for fast dissociating analytes. Using data collection rate 40

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Property	Variable	Value	Flow path	Predip
Solution	<input type="checkbox"/>	Capture solution	2	<input checked="" type="checkbox"/>
Contact time	<input type="checkbox"/>	60 s		
Flow rate	<input type="checkbox"/>	10 µl/min		
Concentration	<input type="checkbox"/>	nM		
Molecular weight	<input type="checkbox"/>	Da		

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

Fig 4.1. Ligand is concentrated on the surface through electrostatic attraction when the pH lies between the isoelectric point of the ligand and the pK_a of the surface. If the pH is too low or too high, ligand will not be concentrated on the surface.

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
Method Builder - 1 Series

1. Method definition | **2. Variables and positioning** | 3. Cycle overview | 4. Plate layout

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2. Send to queueで測定開始

The screenshot shows the 'Method Builder - 1 Series' interface. The '4. Plate layout' step is active, and the 'Send to queue' button is highlighted with a blue box and a circled '3'. The 'Trays' view shows a 12-well plate layout with reagent positions marked. The 'Volume summary' view shows a table of reagent volumes for Rack 1.

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Id	Ligand 1	Ligand 2	50 mM NaOH
4	Immobilization buffer 10 mM acetate pH 4.0 ● 44 µl	Immobilization buffer 10 mM acetate pH 4.0 ● 44 µl	
3	Immobilization buffer 10 mM acetate pH 4.5 ● 44 µl	Immobilization buffer 10 mM acetate pH 4.5 ● 44 µl	
2	Immobilization buffer 10 mM acetate pH 5.0 ● 44 µl	Immobilization buffer 10 mM acetate pH 5.0 ● 44 µl	
1	Immobilization buffer 10 mM acetate pH 5.5 ● 44 µl	Immobilization buffer 10 mM acetate pH 5.5 ● 44 µl	● 282 µl

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Instrument control Methods Runs Action history Preferences

Immobilization x Open + New

Immobilization - 1 Series 1. Method definition 2. Positioning and plate layout Send to queue Compatibility

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View Trays Volume summary ①

Sort positions Ascending Descending ③ Estimated run time 32 min

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
Instrument control Methods **Runs** Action history

LS-044885 immob Anti-B2u low level 4 channels 1 Open

LS-044885 immob Anti-B2u low level 4 channels 1 Results Sensorgrams Run properties

Immobilization result ① ②

Chip CM5

Flow cell	Channel	Procedure	Step	Ligand	Response bound (RU)	Response final (RU)
2	1	Contact time	Amine	Anti-B2u	623.8	440.2
2	2	Contact time	Amine	Anti-B2u	731.3	557.9
2	3	Contact time	Amine	Anti-B2u	745.0	573.0
2	4	Contact time	Amine	Anti-B2u	77.1	-6.0

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



The screenshot shows the 'Interactive run' configuration screen. It includes an 'Activity queue' with 'Interactive run' selected, a 'Checklist' with items like 'Correct chip is docked' and 'Waste bottle is empty', 'Bottles' settings for 'Use buffer position' (A, B, C, D), 'Data collection rate' (1 Hz), 'Use flow cells' (1,2) and 'Reference' (1), 'Concentration units' (nM), and a 'Tray setup' section with a grid of wells (12x4) and a 'Ready to start' button. Red circles with numbers 1-6 are overlaid on the interface to correspond to the table above.

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
① Single cycle kinetics 1
② Regeneration 1

③ Change units

Command type | Concentration unit

Command type	Concentration unit
Analyte	nM
Single cycle kinetics	nM
A-B-A	nM
Dual	nM
Capture	nM
Enhancement	nM
General	nM
Poly	nM

Apply and close | Close

Property | Variable | Value

Solution	<input checked="" type="checkbox"/>	
Contact time	<input type="checkbox"/>	120 s
Dissociation time	<input type="checkbox"/>	600 s
Flow rate	<input type="checkbox"/>	30 μ /min
Molecular weight	<input type="checkbox"/>	Da

Flow path: 1, 2 | Predip: | Concentrations per cycle: 5

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Variables and positioning

Method Builder - 1 Series

1. Method definition | **2. Variables and positioning** | 3. Cycle overview | 4. Plate layout | Send to queue

Settings

Actions for selected step

File | Clipboard

Manage cycles

Add cycle | Insert above | Remove cycles | Remove all | Move up | Move down

Single cycle kinetics 1 -

No	Solution	Control	Concentration 1 (nM)	Concentration 2 (nM)	Concentration 3 (nM)	Concentration 4 (nM)	Concentration 5 (nM)
1	Sample 1	▼	0	0	0	0	0
2	Sample 1	▼	2.4	12	60	300	1500

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
2. Send to queueで測定開始

The screenshot shows the 'Method Builder - 1 Series' interface. The '4. Plate layout' step is selected, and the 'Send to queue' button is highlighted. The 'Bottles' section shows Buffer and Water bottles. The 'Trays' section shows a 96-well plate layout with callout 1. The 'Rack 1' section shows a reagent rack layout with callout 2. The 'Plate 1' table shows sample concentrations and volumes with callout 3.

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
The screenshot shows the software interface with the following elements:

- Header: cytiva Biacore™ Insight Evaluation Software
- Navigation: Create new evaluation, Open existing evaluation, Action history
- Tabs: 1. Select runs (active), 2. Select evaluation method
- Search: Search... (with magnifying glass and close icons)
- Options: Select runs from multiple folders, Show 20 items 1-16
- File Tree (Left):
 - Root
 - Demo Data
 - Biacore 8K Demo data
 - Biacore8K Training data (selected)
 - IE Concentration Demo Files
 - IE Epitope Binning Demo Files
 - IE Intelligent Analysis Demo Files

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The screenshot shows the software interface with the following elements:

- Header: cytiva Biacore™ Insight Evaluation Software
- Navigation: Create new evaluation, Open existing evaluation, Action history
- Tabs: 1. Select runs, 2. Select evaluation method (active)
- Method Selection: User defined, Predefined
- Method List (Left):
 - Empty method
 - Surface preparation
 - Assay development
 - Binding screen
 - Kinetics / affinity
 - Fragment
 - LMW
 - Antibody / general (selected)
- Method List (Right):

Name
Kinetic screen - Evaluation method
Single-cycle kinetics - Evaluation method (selected)
Multi-cycle kinetics - Evaluation method
Multi-cycle affinity - Evaluation method
Parallel kinetics - Evaluation method
2D kinetics - Evaluation method
Single-cycle kinetics using capture - Evaluation method
- Button: Open

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Settings
Data grouping (8)
Serial Parallel/2D
Injection assignment
Kinetics/Affinity mode
Settings apply to selected series
Kinetics Affinity Both
Blank settings
Fit models Perform fit [0]
Initial values
Quality prediction: None

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
3-2. Kinetics/Affinity (K_D, k_a, k_d)測定

Control Softwareの構造～Method definitions

Instrument control | Methods | Runs | Action history | Preferences

Single-cycle kinetics using Biotin CAPture kit | Open | New

Method Builder - 1 Series

1. Method definition | 2. Variables and positioning | 3. Cycle overview | 4. Plate layout | Send to queue

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General settings

Use flow cells: 1, 2
Reference: 1
Change flow cells

Data collection rate: 10 Hz
Running buffer: Buffer
Concentration unit nM
Change units

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2. Send to queueで測定開始

Method Builder - 1 Series

1. Method definition 2. Variables and positioning 3. Cycle overview 4. Plate layout Send to queue

Bottles
Buffer bottle 200 ml Buffer
Water bottle 200 ml Water

View Trays Volume summary


Plate 1

	A	B	C	D	E	F	G	H
10	Sample 1 1500 nM ● 103 μ l							
9	Sample 1 300 nM ● 103 μ l							
8	Sample 1 60 nM ● 103 μ l							
7	Sample 1 12 nM ● 103 μ l							
6	Sample 1 2.4 nM ● 103 μ l							
5	Sample 1 0 nM ● 103 μ l							
4	Sample 1 0 nM ● 103 μ l							
3	Sample 1 0 nM							

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


Settings

Data grouping (8)
Serial Parallel/2D
Injection assignment
Kinetics/Affinity mode
Settings apply to selected series
Kinetics Affinity Both
Blank settings
Fit models Perform fit [0]
Initial values
Quality prediction: None

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