# Biacore 1 series

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

## 3. Insight Control Softwareの起動・ログイン

A screenshot of the Biacore Insight Control Software login interface. The screen displays the Cytiva logo and the software title "Biacore™ Insight Control Software" with version "Version 5.0.16.21762". There are input fields for "User id" (containing "XXXXXX / 12345"), "Domain" (containing "XXXXXX"), and "Password". Below these are options for "Database" (set to "Local Biacore Insight database") and "Selected extensions" (listing "Biacore Intelligent Analysis™", "Concentration & Potency", "Data Integration", "Epitope Binning", and "Extended Screening"). At the bottom, there are "Log in" and "Cancel" buttons.

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
# 1-3. Biacore Insight Control Software

Instrument control画面

The screenshot shows the Biacore Insight Control Software interface. The top navigation bar includes the Cytiva logo, the software name, user information, and a help menu. Below the navigation bar are four main panels, each highlighted with an orange border:

- Activity queue:** Contains a 'Change solutions' task with a 6-minute progress bar and an 'Interactive run' button. An 'Add activity' button is at the bottom.
- System setup:** Lists tasks such as 'Set flow cell temperature', 'Set sample compartment temperature', 'Change chip', 'Change solutions', 'Immobilization checkpoint', and 'Wait'.
- Maintenance:** Lists tasks such as 'Desorb', 'Desorb and sanitize', 'System check', 'Normalize', and 'Shutdown'.
- Start:** Contains 'Interactive run' and 'Method' buttons.

Below the panels, there are four descriptive labels in green text:

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The bottom status bar shows system information: 'Running Change solutions...', 'Current buffer Buffer', 'CAP: 3/3/2023 6:09:55 PM', 'chip information', 'Flow cell 25.0 °C (25)', 'Sample compartment 25.0 °C (25)', 'Illumination on/off', and 'Hotel door is closed'.

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
The screenshot shows the Biacore Insight Control Software interface. The 'Methods' tab is selected in the top navigation bar. Below the navigation bar, there are 'Open' and 'New' buttons. The 'New' button is highlighted with a red box. The main content area is divided into two sections: 'Biacore 8 series' and 'Biacore 1 series'. Under 'Biacore 1 series', the 'Kinetics / affinity' folder is expanded, showing a list of method templates. A table displays the details of these templates.

Name	Description
Kinetic screen	Single concentration injections for rapid estimation of kinetic properties. Includes...
Single-cycle kinetics	Single-cycle kinetics for protein samples.
Multi-cycle kinetics/affinity	Multi-cycle kinetics/affinity for protein samples. Includes regeneration.
Single-cycle kinetics using capture	Single-cycle kinetics for protein samples using captured ligand.
Multi-cycle kinetics/affinity using capture	Multi-cycle kinetics/affinity for protein samples, using captured ligand. Includes...
Fast kinetics 10 Hz	Multi-cycle kinetics for fast dissociating analytes. Using data collection rate 10
Fast kinetics 40 Hz	Multi-cycle kinetics for fast dissociating analytes. Using data collection rate 40

On the right side of the interface, there is an 'Item details' panel with fields for 'Name' and 'Description'. An 'Open' button is located at the bottom right of the interface.

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General settings  
Use flow cells 1, 2  
Reference 1  
Change flow cells

Data collection rate  
1 Hz

Running buffer  
Buffer

Concentration unit varies  
Change units

Name Samples  
Purpose Analysis

Startup Samples Solvent correction Negative control Positive control Add step

④ ⑤ ⑥

Analyte 1 Wash 1 Carry over control 1 Add command

① ② ③

Property Variable Value Type  
Solution  Flow path 1, 2  
Contact time  60 s  
Dissociation time  60 s  
Flow rate  30 µl/min  
Concentration  µM  
Molecular weight  Comment

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Variables and positioning

Method Builder - 1 Series

1. Method definition 2. Variables and positioning 3. Cycle overview 4. Plate layout

Settings

Actions for selected step

File

Clipboard

Manage cycles

Add cycle

Insert above

Remove cycles

Remove all

Move up

Move down

25 °C

Startup

Samples

Solvent correction

Negative control

Positive control

No	Analyte 1 - Solution	Control	Molecular Weight (Da)
1	Sample 1		324
2	Sample 2		353
3	Sample 3		567
4	Sample 4		346
5	Sample 5		75
6	Sample 6		678
7	Sample 7		1063
8	Sample 8		765
9	Sample 9		479
10	Sample 10		598
11	Sample 11		789
12	Sample 12		546
13	Sample 13		444
14	Sample 14		243
15	Sample 15		573
16	Sample 16		732
17	Sample 17		729

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
2. Send to queueで測定開始

The screenshot shows the 'Plate layout' step in the software. The 'Trays' tab is active, showing a 96-well plate layout with samples and reagents. The 'Volume summary' tab is also visible. A table below explains the interface elements:

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
1. Select runs 2. Select evaluation method

Search...

Select runs from multiple folders Show 20 items 1-12 of 12

- Root
  - Demo Data
    - 1 series
    - Biacore 8K Demo data
      - v2**
      - v4
    - Biacore8K Training data
    - IE Concentration Demo Files
    - IE Epitope Binning Demo Files

Name
Regeneration scouting 6/21/2016 12:19:11 PM
Clean screen (large) 2016-06-03
Affinity Screen Analogues 2016-05-04
Affinity screen 2016-04-27 17:39:29
BLS (small) 2016-04-24
Clean Screen (small) 2016-04-24
<b>BLS (large) 20160316</b>
Fragment screen GX

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1. Select runs 2. Select evaluation method

User defined Predefined

- Empty method
- Surface preparation
- Assay development
- Binding screen
  - Fragment
  - LMW**
  - Antibody / general
- Kinetics / affinity

Name
LMW screen - Evaluation method
LMW screen using capture - Evaluation method

Open

# 3-2. 低分子のBinding screen解析

Solvent correction

Biacore Insight Evaluation Software - BLS (small) 2016-04-24

cytiva Biacore™ Insight Evaluation Software User HELP

Create new evaluation Open existing evaluation Evaluation - BLS (small) 2016-04-24 Action history

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# 3-2. 低分子のBinding screen解析

Evaluation - Binding 画面 (解析結果)

Biocore Insight Evaluation Software - BLS (large) 20160316

cytiva Biacore™ Insight Evaluation Software User HELP

Create new evaluation Open existing evaluation Evaluation - BLS (large) 20160316 Action history

QC - Sensorgram QC - Baseline QC - Binding to reference Evaluation - Binding Home

Settings

- Data grouping: All together (1)
- Injection assignment
- Axis settings
- Adjustments: Applied
- Boundaries: Control based
- Binder prediction: None
- Color by

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Cycle	Channel	Sensorgram type	Analysis step purpose	Analysis step name	Excluded	Curve markers	Analyte 1 Solution	Analyte binding late_1 Adjusted relative	Analyte 1 Control type	Cut-off
81	8	Corrected	Analysis	Samples and Controls			Filibuvir	101.3	Positive control	Above
82	8	Corrected	Analysis	Samples and Controls			DMSO	22.56	Negative control	Below
83	8	Corrected	Analysis	Samples and Controls			DMSO	23.99	Negative control	Below
84	8	Corrected	Analysis	Samples and Controls			DMSO	26.04	Negative control	Below
85	8	Corrected	Analysis	Samples and Controls		Solvent corr. out of range	FRAG361		Not a control	
86	8	Corrected	Analysis	Samples and Controls			FRAG368	154.3	Not a control	Above
87	8	Corrected	Analysis	Samples and Controls			FRAG375	948.2	Not a control	Above
88	8	Corrected	Analysis	Samples and Controls			FRAG389	76.24	Not a control	Above

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
# 4-3. Kinetics/Affinity ( $K_D$ , $k_a$ , $k_d$ )測定

Instrument control Methods Runs Action history Preferences

LMW single-cycle kinetics using Biotin CAPture kit Open New

Method Builder - 1 Series

1. Method definition 2. Variables and positioning 3. Cycle overview 4. Plate layout Send to queue

25 °C Conditioning Startup Analysis Solvent correction First/Every 48 cycles/Last Add step

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⑦

Command type Concentration unit

Command type	Concentration unit
Analyte	nM
Single cycle kinetics	nM
A-B-A	nM
Dual	nM
Capture	nM
Enhancement	nM
General	nM
Poly	nM

Apply and close Close

Name: Analysis Purpose: Analysis

Flow cell 1 Ref General 1 Capture 1 Single cycle kinetics 1 Regeneration 1 Wash 1 Wash 2

Property	Variable	Value
Solution	<input type="checkbox"/>	Biotin CAPture
Contact time	<input type="checkbox"/>	300 s
Dissociation time	<input type="checkbox"/>	0 s
Flow rate	<input type="checkbox"/>	2 $\mu$ l/r
Concentration	<input type="checkbox"/>	nM
Molecular weight	<input type="checkbox"/>	Da

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# 4-3. Kinetics/Affinity ( $K_D$ , $k_a$ , $k_d$ )測定

Variables and positioning

The screenshot shows the 'Method Builder - 1 Series' interface. The '2. Variables and positioning' tab is active. It contains a table for 'Single cycle kinetics 1' and a well plate layout for 'Plate 1' and 'Rack 1'. The table has columns for 'No', 'Solution', 'Control', and five concentration columns. The well plate layout shows a 12x8 grid for Plate 1 and a 12x4 grid for Rack 1. A 'Positioning settings' sidebar is visible on the right.

No	Solution	Control	Concentration 1 (nM)	Concentration 2 (nM)	Concentration 3 (nM)	Concentration 4 (nM)	Concentration 5 (nM)
1	Sample 1	▼	0	0	0	0	0
2	Sample 1	▼	2.4	12	60	300	1500

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2. Send to queueで測定開始

Method Builder - 1 Series

1. Method definition 2. Variables and positioning 3. Cycle overview 4. Plate layout Send to queue

Bottles

Buffer bottle 200 ml Buffer

Water bottle 200 ml Water

View Trays Volume summary


Plate 1

	A	B	C	D	E	F	G	H
10	Sample 1 1500 nM ● 103 $\mu$ l							
9	Sample 1 300 nM ● 103 $\mu$ l							
8	Sample 1 60 nM ● 103 $\mu$ l							
7	Sample 1 12 nM ● 103 $\mu$ l							
6	Sample 1 2.4 nM ● 103 $\mu$ l							
5	Sample 1 0 nM ● 103 $\mu$ l							
4	Sample 1 0 nM ● 103 $\mu$ l							
3	Sample 1 0 nM ● 103 $\mu$ l							

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
1. Select runs 2. Select evaluation method

Search...

Select runs from multiple folders Show 20 items 1-16

- Root
  - Demo Data
    - Biacore 8K Demo data
      - Biacore8K Training data**
      - IE Concentration Demo Files
      - IE Epitope Binning Demo Files
      - IE Intelligent Analysis Demo Files

Name
System check 6/9/2017 10:29:35 AM
SCK demo 11/30/2016 1:23:47 PM
Fragment affinity screen 11/2/2016 3:33:31 PM
Fragment clean screen 11/2/2016 12:31:04 PM
Immobilization Clean Screen 11/2/2016 11:09:57 AM
<b>Single-cycle kinetics 011116 11/1/2016 3:40:51 PM</b>

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1. Select runs 2. Select evaluation method

User defined **Predefined**

- Empty method
- Surface preparation
- Assay development
- Binding screen
- Kinetics / affinity**
  - Fragment
  - LMW**
  - Antibody / general

Name
LMW single-cycle kinetics - Evaluation method
LMW multi-cycle kinetics - Evaluation method
LMW multi-cycle affinity - Evaluation method
<b>LMW single-cycle kinetics using capture - Evaluation method</b>
LMW multi-cycle kinetics using capture - Evaluation method
LMW multi-cycle affinity using capture - Evaluation method

Open

# 4-4. Kinetics/Affinity ( $K_D$ , $k_a$ , $k_d$ )解析

Solvent correction

Biacore Insight Evaluation Software - 2\_CA\_drug

cytiva Biacore™ Insight Evaluation Software User ? HELP

Create new evaluation Open existing evaluation Evaluation - 2\_CA\_drug Action history

Solvent correction

Included	#Excluded	Cycle	Flow cell	Chi <sup>2</sup> (RU <sup>2</sup> )	Y0 (RU)	Immobilization level (RU)
<input checked="" type="checkbox"/>	0	5	4-3	0.427	0.9	
<input checked="" type="checkbox"/>	0	18	4-3	0.672	2.2	
<input checked="" type="checkbox"/>	0	31	4-3	0.929	2.5	
<input checked="" type="checkbox"/>	0	44	4-3	0.397	2.6	
<input checked="" type="checkbox"/>	0	57	4-3	0.308	3.1	
<input checked="" type="checkbox"/>	0	70	4-3	0.284	3.4	

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**Settings**  
Data grouping (8)  
Serial Parallel/2D  
Injection assignment  
Kinetics/Affinity mode  
Settings apply to selected series  
Kinetics Affinity Both  
Blank settings  
Fit models Perform fit [0]  
Initial values  
Quality prediction: None

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