
Biacore 1 series

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Version 1.2
2023/09

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


3. Insight Control Softwareの起動・ログイン

A screenshot of the Biacore Insight Control Software login interface. The interface is white with a 'cytiva' logo at the top. Below the logo, it says 'Biacore™ Insight Control Software' and 'Version 5.0.16.21762'. There are input fields for 'User id' (containing 'XXXXXX / 12345'), 'Domain: XXXXX', and 'Password'. Below these are sections for 'Database' (set to 'Local Biacore Insight database') and 'Selected extensions' (listing 'Biacore Intelligent Analysis™', 'Concentration & Potency', 'Data Integration', 'Epitope Binning', and 'Extended Screening'). At the bottom, there are 'Log in' and 'Cancel' buttons.

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
3. Biacore Insight Control Software

Instrument Control画面

The screenshot shows the Biacore Insight Control Software interface. The top navigation bar includes the Cytiva logo, the software name, user information, and a help menu. Below this is a secondary navigation bar with tabs for 'Instrument control', 'Methods', 'Runs', and 'Action history'. The main content area is divided into four primary sections: 'Activity queue', 'System setup', 'Maintenance', and 'Start'. Each section contains a list of actions with icons and a progress indicator. The 'Activity queue' section has a progress bar showing 6 minutes. The 'System setup' section includes actions like 'Set flow cell temperature', 'Set sample compartment temperature', 'Change chip', 'Change solutions', 'Immobilization checkpoint', and 'Wait'. The 'Maintenance' section includes 'Desorb', 'Desorb and sanitize', 'System check', 'Normalize', and 'Shutdown'. The 'Start' section includes 'Interactive run' and 'Method'. A status bar at the bottom provides real-time information such as 'Running Change solutions...', 'Current buffer', 'chip information', 'Flow cell' temperature (25.0 °C), 'Sample compartment' temperature (25.0 °C), 'Illumination' status (on/off), and 'Hotel door is closed'.

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


The screenshot shows the 'Biacore™ Insight Control Software' interface. The top navigation bar includes 'Instrument control', 'Methods', 'Runs', and 'Action history'. The main content area is divided into several sections: 'Activity queue' with an 'Interactive run' entry marked 'Action required'; 'Checklist' with items like 'Correct chip is docked' and 'Waste bottle is empty'; 'Bottles' with instructions for buffer and reagent tubes; 'Data collection rate' set to 1 Hz; 'Use channels' with checkboxes for channels 1-8 and 'Check All'; and 'Concentration units' set to nM. On the right, the 'Tray setup' section shows a 96-well plate grid and a 'Type' dropdown set to 'None'. A 'Ready to start' button is visible at the bottom.

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Property	Variable	Value	Flow path	Predip
Solution	<input type="checkbox"/>	Capture solution	2	<input checked="" type="checkbox"/>
Contact time	<input type="checkbox"/>	60 s		
Flow rate	<input type="checkbox"/>	10 μ /min		
Concentration	<input type="checkbox"/>			
Molecular weight	<input type="checkbox"/>			

8. Kinetics/Affinity (K_D , k_a , k_d)測定

Variables and positioning

Single-cycle kinetics using capture

Method Builder - 1 Series

1. Method definition 2. Variables and positioning 3. Cycle overview 4. Plate layout

Settings

Actions for selected step

File

Clipboard

Manage cycles

Add cycle

Insert above

Remove cycles

Remove all

Move up

Move down

No	Solution	Control	Concentration 1 (nM)	Concentration 2 (nM)	Concentration 3 (nM)	Concentration 4 (nM)	Concentration 5 (nM)
1	Sample 1		0	0	0	0	0
2	Sample 1		2.4	12	60	300	1500

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8. Kinetics/Affinity (K_D, k_a, k_d)測定

Positioning settings

Single-cycle kinetics using capture

Method Builder - 1 Series

1. Method definition 2. Variables and positioning 3. Cycle overview 4. Plate layout

Settings

Actions for selected step

Import from

File

Clipboard

Manage cycles

Add cycle

Insert above

Remove cycles

Remove all

Move up

Move down

Positioning settings

Groups in priority order	Pooling	Plate Rack	Vial size	Priority
Single cycle kinetics 1 Step Analysis	Yes No		Small Large	
Capture 1 Step Startup	Yes No		Small Large	
Capture 1 (Predip) Step Startup				
Capture 1 Step Analysis				
Capture 1 (Predip) Step Analysis				
Analyte 1 Step Startup				
Regeneration 1				

Plate and Rack

Plate positioning direction

Horizontal Vertical

Fill both plate and rack before starting the next tray

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

9. Biacore Insight Evaluation Software

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The screenshot shows the 'Select runs' step in the software. The left sidebar shows a file tree with 'Biacore8K Training data' selected. The main area displays a list of runs with columns for Name and Date. The selected run is 'Single-cycle kinetics 011116 11/1/2016 3:40:51 PM'. A blue button labeled 'Select evaluation method' is positioned at the bottom center.

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The screenshot shows the 'Select evaluation method' step. The '2. Select evaluation method' tab is active. The 'Predefined' section is expanded, and 'Single-cycle kinetics - Evaluation method' is selected. A blue button labeled 'Open' is positioned at the bottom right.

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