
# Biacore 8 series

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Version 1.2  
2023/09

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# 3. Biacore Insight Control Software

Instrument Control画面

The screenshot shows the Biacore Insight Control Software interface. The top navigation bar includes the Cytiva logo, the software name, user information, and a help menu. Below this is a secondary navigation bar with tabs for 'Instrument control', 'Methods', 'Runs', and 'Action history'. The main workspace is divided into four primary sections: 'Activity queue', 'System setup', 'Maintenance', and 'Start'. Each section contains a list of actions with corresponding icons. The 'Activity queue' section includes 'Change solutions' and 'Interactive run'. 'System setup' includes 'Set flow cell temperature', 'Set sample compartment temperature', 'Change chip', 'Change solutions', 'Immobilization checkpoint', and 'Wait'. 'Maintenance' includes 'Desorb', 'Desorb and sanitize', 'System check', 'Normalize', and 'Shutdown'. The 'Start' section includes 'Interactive run' and 'Method'. A status bar at the bottom provides real-time information such as 'Running Change solutions...', 'Current buffer', 'chip information', 'Flow cell' temperature (25.0 °C), 'Sample compartment' temperature (25.0 °C), 'Illumination' status (on/off), and 'Hotel door is closed'.

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


The screenshot shows the 'Biacore™ Insight Control Software' interface. The top navigation bar includes 'Instrument control', 'Methods', 'Runs', and 'Action history'. The main area is divided into several sections: 'Activity queue' with an 'Interactive run' entry marked 'Action required'; 'Checklist' with items like 'Correct chip is docked' and 'Waste bottle is empty'; 'Bottles' with instructions for buffer and reagent tubes; 'Data collection rate' set to 1 Hz; 'Use channels' with checkboxes for channels 1-8 and 'Check All'; and 'Concentration units' set to nM. On the right, the 'Tray setup' section shows a 96-well plate grid and a 'Type' dropdown set to 'None'. A blue 'Ready to start' button is at the bottom.This screenshot shows a close-up of a dropdown menu. The menu is open, showing two options: 'Flow cell 1' and 'Flow cell 2'. A dark grey button with the text 'Add command' and a downward-pointing chevron is positioned over the menu.

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# 8. Kinetics/Affinity ( $K_D, k_a, k_d$ )測定

Control Softwareの構造～Method definitions

The screenshot shows the Biacore Insight Control Software Method Builder interface. The top navigation bar includes 'Instrument control', 'Methods', 'Runs', and 'Action history'. The main workspace is titled 'Method Builder - 8 Series' and contains four tabs: '1. Method definition', '2. Variables and positioning', '3. Cycle overview', and '4. Plate layout'. The 'Method definition' tab is active, showing a flow diagram with steps: 'Startup', 'Analysis', 'Single cycle kinetics 1', and 'Regeneration 1'. The 'Analysis' step is highlighted with a green arrow and a text box: 'Method definitionにおけるフロー Method definitionステップ内の各Assay stepのタブ'. Below the flow diagram, the 'Analysis' step is expanded to show a table of properties and values, and an 'Add command' button is highlighted with an orange arrow and a text box: 'Analysisにおける各コマンド 各Assay step (この場合"Analysis") 内の添加条件などの設定'. On the left, a 'General settings' panel is visible. Annotations in blue and orange point to specific parts of the interface.

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The screenshot displays the Biacore Insight Control Software interface. The main window is titled "Single-cycle kinetics using capture" and shows a "Method Builder - 8 Series" configuration. The "Variables and positioning" step is selected in the top navigation bar. The interface is divided into several sections:

- Settings:** Includes "Use channels" (1-8) and "Actions for selected step" (File, Clipboard, Epitope binning, Add cycle, Insert above, Remove cycles, Remove all, Move up, Move down).
- Startup and Analysis:** Two tabs are visible. The "Analysis" tab is active, showing a table for "Single cycle kinetics 1 -".
- Positioning:** Shows four plate layouts (Plate 1, Plate 2, Plate 3, Plate 4) with well positions (A-H, 1-12) and a "Positioning settings" panel.

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No	Solution	Control	Concentration 1 (nM)	Concentration 2 (nM)	Concentration 3 (nM)	Concentration 4 (nM)	Concen
1	Sample 1	▼	0	0	0	0	0
	Sample 2	▼	0	0	0	0	0
	Sample 3	▼	0	0	0	0	0
	Sample 4	▼	0	0	0	0	0
	Sample 5	▼	0	0	0	0	0
	Sample 6	▼	2.4	12	60	300	
	Sample 7	▼	2.4	12	60	300	
	Sample 8	▼	2.4	12	60	300	

# 8. Kinetics/Affinity ( $K_D$ , $k_a$ , $k_d$ )測定

Positioning settings

The screenshot displays the Biacore Insight Control Software interface. The top navigation bar includes the Cytiva logo, software name, user name (Gen Takata), and instrument ID (Biacore 8K 1001162636). The main menu shows 'Instrument control', 'Methods', 'Runs', and 'Action history'. The current method is 'Single-cycle kinetics using capture'. The 'Method Builder - 8 Series' window is open, showing four tabs: '1. Method definition', '2. Variables and positioning', '3. Cycle overview', and '4. Plate layout'. The '2. Variables and positioning' tab is active, showing a list of samples and two plate layouts (Plate 1 and Plate 2). The 'Positioning settings' panel is open, showing a list of steps with 'Pooling', 'Plate', and 'Priority' options. The 'Pooling' option is checked for 'Single cycle kinetics 1 Step Analysis' and 'Capture 1 Step Startup'. The 'Plate' option is set to 'LR' for 'Capture 1 Step Startup'. The 'Priority' option is set to 'LR' for 'Capture 1 Step Startup'. The 'Tray' section shows 'Number of trays' set to 'Auto populate'. The 'Series' section shows 'Keep series together' checked. The 'Reposition' section shows 'Autopopulate position after manual reposition' checked. The 'Actions' section is empty.

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

# 9. Biacore Insight Evaluation Software

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The screenshot shows the '1. Select runs' step in the software. The left sidebar displays a file tree with 'Biacore8K Training data' selected. The main area shows a list of runs, with 'Single-cycle kinetics 011116 11/1/2016 3:40:51 PM' highlighted. A blue button labeled 'Select evaluation method' is positioned at the bottom center.

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The screenshot shows the '2. Select evaluation method' step. The 'Predefined' tab is active, and 'Single-cycle kinetics - Evaluation method' is selected in the list. A blue button labeled 'Open' is positioned at the bottom right.

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