
# Biacore T200/S200

## Insight Evaluation Software 操作手順書

Ver.202108

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
The screenshot displays the Biacore Insight Evaluation Software interface. The main window is titled "Biacore™ Insight Evaluation Software" and has two tabs: "Create new evaluation" and "Open existing evaluation". The "Open existing evaluation" tab is active, showing a workflow with two steps: "1. Select runs" and "2. Select evaluation method".


Step 2, "2. Select evaluation method", is currently selected. It features a sidebar with a tree view of methods. The "Predefined" section is expanded, and the "LMW" method is highlighted. A list of predefined methods is shown in the center, including "LMW screen - Evaluation method" and "LMW screen using capture - Evaluation method". The "LMW screen - Evaluation method" is selected.

On the right side of the interface, there are two panels. The top panel, "Selected runs in order", lists a run: "BLS (large) 20160316" with its path and description. The bottom panel, "Selected evaluation method", shows the selected "LMW screen - Evaluation method" with its description: "Evaluation of single concentration LMW screen. Identifies hits based on binding level." At the bottom right, there is an "Open" button.

Numbered callouts (② through ⑥) are placed on the screenshot to correspond with the steps in the table above: ② points to the "1. Select runs" tab, ③ to the "2. Select evaluation method" tab, ④ to the "Predefined" section in the sidebar, ⑤ to the "LMW" method in the sidebar, ⑥ to the "Open" button, and an additional ⑤ points to the selected "LMW screen - Evaluation method" in the center panel.

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




The screenshot shows the Biacore Insight Evaluation Software interface. The main window title is "Biacore™ Insight Evaluation Software". At the top, there are buttons for "Create new evaluation" and "Open existing evaluation". Below this, there are two tabs: "1. Select runs" and "2. Select evaluation method". The "2. Select evaluation method" tab is active, showing a list of predefined methods under the "Predefined" category. The list includes various kinetic methods, with "LMW single-cycle kinetics using capture - Evaluation method" selected. To the right, there is a "Selected runs in order" panel showing a single run: "Single-cycle kinetics 011116 11/1/2016 3:40:51 PM". Below that, the "Selected evaluation method" panel shows the details for the selected method: "LMW single-cycle kinetics using capture - Evaluation method" with a description: "Evaluation of single-cycle kinetics for LMW samples binding to captured ligand, to calculate kinetic rate constants." At the bottom right, there is an "Open" button. Numbered annotations 1 through 6 are placed on the interface to correspond to the steps in the table above.

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The screenshot displays the Biacore Insight Evaluation Software interface. The main window is titled "Biacore™ Insight Evaluation Software" and shows a "Solvent correction" dialog box. The dialog box contains a table with columns: Included, #Excluded, Cycle, Channel, CH<sup>2</sup> (RU<sup>2</sup>), YO (RU), and Immobilization level (RU). The table lists 24 rows of data, all with "Included" checked. To the right of the table is a graph showing "Response (RU) vs Response Reference (RU)" with multiple curves. Below the graph is another graph showing "Response (RU) vs Time (s)" with a step-like signal. At the bottom left, there is a button labeled "Apply and close" with a circled 2 next to it, and a "Cancel" button.

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

The screenshot shows a software interface with a table of analysis steps. The table has columns for Cycle, Channel, Analysis step purpose, Analysis step name, Analyte Solution, and Analyte Concentration. A dropdown menu is open over the 'Analyte Solution' column of the 4th row, showing options: Not a control, Positive control, Negative control, and Blank control.

Run number	Concentration unit
1	nM

Cycle	Channel	Analysis step purpose	Analysis step name	Analyte Solution	Analyte Concentration
1	1	Startup	Startup	Buffer	
2	1	Startup	Startup	Buffer	
3	1	Startup	Startup	Buffer	
4	1	Analysis	Sample	b2m	
5	1	Analysis	Sample	b2m	
6	1	Analysis	Sample	b2m	
7	1	Analysis	Sample	b2m	
8	1	Analysis	Sample	b2m	
9	1	Analysis	Sample	b2m	

- Not a control
- Positive control
- Negative control
- Blank control

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The screenshot displays the 'Run 1' configuration page. At the top, there is a 'Home' button and a 'Run 1' tab, with an arrow pointing to the tab labeled ④. Below this, the page is split into two columns: 'Chip information' and 'Run information'. The 'Chip information' column lists details like 'Chip type: CM5', 'Chip id: 210531-1147:2623147', 'Lot number', and 'First dock date: 5/31/2021 11:51:12 AM'. The 'Run information' column lists 'Name: 210531\_test\_run01', 'Path: Root\Guest\test\_run01', 'Method', and 'Instrument type: Biacore S200'. A 'More info' button is located to the right of the 'Run information' section. Below these sections are two 'Flow cell' configuration areas. 'Flow cell 1' has a 'Channel 1' checkbox checked and a 'Sensorgram type' dropdown menu with 'Active' and 'Reference' options; an arrow points to the 'Reference' option labeled ⑤. 'Flow cell 2' is partially visible below it. At the bottom of the page, there are 'Apply and close' and 'Cancel' buttons, with an arrow pointing to the 'Cancel' button labeled ⑤.

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