

仕様				
タンパク質	分子量 (Da)	由来	参考文献	タンパク質/バイアル (μg)
HMW-SDS (53,000-220,000)				
Myosin	220,000	rabbit muscle	1	25
α-2-Macroglobulin	170,000	bovine plasma	2	100
β-Galactosidase	116,000	<i>E. coli</i>	3	16
Transferrin	76,000	human	4	17
Glutamate dehydrogenase	53,000	bovine liver	5	18
LMW (14,400-97,000)				
Phosphorylaseb	97,000	rabbit muscle	6	67
Albumin	66,000	bovine serum	7	83
Ovalbumin	45,000	chicken egg white	8	147
Carbonic anhydrase	30,000	bovine erythrocyte	9	83
Trypsin inhibitor	20,100	soybean	10	80
α-Lactalbumin	14,400	bovine milk	11	116
HMW Native (66,000-669,000)				
Thyroglobulin	669,000	porcine thyroide	12	76
Ferritin	440,000	equine spleen	13	50
Catalase	232,000	bovine liver	14	36
Lactate dehydrogenase	140,000	bovine heart	15	48
Albumin	66,000	bovine serum	7	40