



















## ご注文情報

Illustration	製品	材質 *1	Applications	包装	コード番号
<b>TUBING CONNECTORS</b>					
	Tubing connector, M6 male for o.d. 1.6 and 1.8 mm tubing	PP	FPLC, for 1.6 and 1.8 mm o.d. tubing	5 個	19-7476-01
	Tubing connector, M6 male for o.d. 1/16" tubing (i.d. 0.5 mm)	PP	FPLC, for 1/16" o.d. tubing	5 個	18-3886-01
	Tubing End fitting required for SRTC *2	PEEK	有機溶媒耐性	1 個	19-2148-01
	Tubing connector, M6 male (2.7 mm o.d.)	PP	Pump P-6000, for 2.7 mm o.d. tubing	5 個	18-4652-01
	Tubing connector, P-1	PP	for 1.8 mm o.d. tubing	6 個	19-2150-01
	Tubing connectors flangeless / M6 female	PP	for 1.8 mm o.d. tubing	2 個	18-1003-68
	Tubing connectors flangeless / M6 male	PP	for 1.8 mm o.d. tubing	5 個	19-1993-50
	Tubing connector TC-2	PP	for 1.8 mm o.d. tubing	2 個	19-0045-01
	Tubing connector TC-3	PP	for 1.8 mm o.d. tubing	1 個	19-0047-01
	Fingertight connector, 1/16" male	PEEK	for 1/16" o.d. tubing	10 個	18-1112-55
	Connector 1/8" male Ferrules, 1/8" (黄)	PEEK PEEK	for 1/8" o.d. tubing	10 個 10 個	18-1121-17 18-1121-18
	Connector 3/16" male Ferrules, 3/16" (青)	PEEK PEEK	for 3/16" o.d. tubing	10 個 10 個	18-1112-49 18-1112-48
	Connector 1/16" male Ferrules, 1/16" (青)	PEEK PEEK	for 1/16" o.d. tubing	10 個 10 個	18-1127-07 18-1127-06
	SRTC-3, M6 female (1.5 mm i.d.)	PEEK		2 個	19-2144-01
	SRTC-2, M6 female (0.5 mm i.d.)	PEEK	2 HR columns in series	1 個	18-3856-01
	SRTC-2, M6 female (1.5 mm i.d.)	PEEK	2 XK columns in series	5 個	19-2143-01
	HiTrap/HiPrep 1/16" male connector for AKTAdesign	PP	to connect top of HiTrap and HiPrep columns to AKTAdesign	8 個	28-4010-81
	Fingertight connector, 1/16" male, narrow head	PEEK	for 1/16" o.d. tubing narrow head	10 個	18-1172-63

\* 1 PP = Polypropylene, PEEK = Polyetheretherketone, SS = Stainless Steel, PA = polyamide, PFR = Perfluororubber

\* 2 Washer 25 個、O-ring 15 個、End fitting 10 個が入っています。

## ご注文情報

Illustration	製品	材質 *1	Applications	包装	コード番号
<b>UNIONS</b>					
	Union M6 female / 1/16" male (Valco compatible, PEEK)	PEEK	to connect columns with M6 connectors to ÄKTAdesign	5 個	18-3858-01
	Union 1/16" female / M6 male	PEEK	to connect columns with 1/16" fittings to FPLC System	6 個	18-1112-57
	Union M6 female / 1/16" male (ÄKTA compatible)	PEEK	to connect columns with M6 connectors to ÄKTAdesign	8 個	18-1112-58
	Union 5/16" female / M6 male	PEEK		3 個	18-1127-76
	Union 1/16" female / 1/16" female (Titanium)	Titanium		1 個	18-3855-01
	Union 1/16" female / 1/16" female (PEEK) *2	PEEK		5 個	11-0003-39
	Union 1/16" male / 1/16" male (0.25 mm i.d.)	PEEK		2 個	18-1120-92
	Union 1/16" male / 1/16" male (0.5 mm i.d.)	PEEK		5 個	28-9543-26
	Union M6 female / 1/16" female (Valco compatible, titanium)	Titanium	HR column in HPLC system	1 個	18-3859-01
	Union M6 female / 1/16" female (Valco compatible, PEEK) *3	PEEK	HR column in HPLC system	2 個	18-1123-94
	Union 5/16" female / 1/16" male	PEEK	to connect columns with 5/16" connectors to ÄKTAdesign	8 個	18-1142-08
	Union M6 female / 5/16" male	PEEK		2 個	18-1169-16
	Union M6 female / 5/16" female	PEEK		2 個	18-1169-17
	Union 5/16" female / 5/16" female	PEEK		2 個	18-1173-51
<b>STOP PLUGS</b>					
	Domed Nut M6	PP	Sealing ends of HR and XK columns	4 個	18-2450-01
	Stop plug 1/16" female	PP	Sealing ends of HiTrap columns	5 個	11-0004-64
	Stop plug M6	PP	Pump P-500 pressure test and to seal unused valve and mixer ports	2 個	19-5170-01
	Fingertight Stop plug, 1/16" male, narrow head	PEEK	Stop plug for use within ÄKTAdesign	5 個	11-0003-55
	Stop plug, 5/16" male	PEEK	Stop plug for use within ÄKTAdesign	5 個	18-1112-50
<b>VALVES</b>					
	SRV-1 (有機溶媒耐性) (ストップコック)	ETFE		1 個	19-2145-01

\* 1 PP = Polypropylene, PEEK = Polyetheretherketone, SS = Stainless Steel, ETFE = Ethylenetetrafluoroethylene, PA = polyamide, PFR = Perfluororubber

\* 2 valco 形式の HiLoad カラム、XK カラムを ÄKTAdesign に接続し、配管延長が必要な場合に使用します。

\* 3 M6 規格の HiLoad カラム、XK カラムを ÄKTAdesign に接続し、配管延長が必要な場合に使用します。

## ご注文情報

Illustration	製品	材質*	Applications	包装	コード番号
<b>LUER CONNECTORS</b>					
	Union adaptors Luerlock female / M6 male	PP	Injection valves V-7, IV-7, MV-7 and IMV-7	2 個	18-1027-62
	Union, Luer female / M6 female	PP		2 個	18-1027-12
	1/16" male/Luer female	PEEK	Injection valve INV-907	2 個	18-1112-51
<b>TCコネクター</b>					
	Connector TC25 / M6 female	PP		2 個	18-1031-09
	Connector TC25 / 5/16" female	PEEK		2 個	18-1169-22
	Connector TC25 / 5/16" male	PEEK		2 個	18-1169-23
	Connector TC25 / 1/4" threaded	PP		2 個	18-0251-98
	Clamp TC25	PA		4 個	18-1169-18
	TC gasket 25/4 mm, PFR	PFR		2 個	18-1169-24
	TC gasket 25/6.5 mm, PFR	PFR		4 個	18-1169-25

\* PP = Polypropylene, PEEK = Polyetheretherketone, SS = Stainless Steel, ETFE = Ethylenetetrafluoroethylene, PA = polyamide, PFR = Perfluororubber