

2-D Quant Kit と併用できる試薬の例

試薬	濃度
SDS	2 % (w/v)
CHAPS	4 % (w/v)
Triton X-100	1 % (w/v)
Pharmalyte pH 3-10	2 % (v/v)
IPG Buffer pH 3-10 NL	2 % (v/v)
Tris	50 mM
EDTA	10 mM
DTT	1 % (65 mM)
2-Mercaptoethanol	2 % (v/v)
Urea	8 M
Thiourea	2 M
Glycerol	30 % (w/v)