

仕様

含まれるタンパク質	分子量 (Mr)	由来	内容量 (およその μg 数)
Pre-Stained Indicator Proteins			
Phosphorylase b	100,000	Rabbit Muscle	50
Myoglobin	16,000	Horse Muscle	50
Lysozyme	15,000	Chicken Egg White	50
Recombinant Tagged Proteins			
Recombinant protein	150,000		1
Recombinant protein	100,000		1
Recombinant protein	75,000		1
Recombinant protein	50,000		1
Recombinant protein	35,000		1
Recombinant protein	25,000		1
Recombinant protein	15,000		1